

# Instruction for Use

**Product Name:** Disposable sampling swab

**Model:** GW-1237NP、GW1237AR、GW1237R

**Product introduction:** Disposable sampling swab is made of tip of flocking nylon fiber, rayon, polyester fiber, PP rod, PS rod and ABS rod. It is mainly suitable for sampling.

**Usage:** The product is sterilized by ethylene oxide, used for clinical sampling.

## Notes:

1. The product is ETO sterilized. Shelf life is 3 years
2. The product is for single use only.
3. Do not use it if expired or package is damaged.

## Instructions:

- 1) Firstly, check whether the outer packaging of a single package is damaged or leaking. If there is such a situation, it is prohibited to use it.
- 2) When there is no abnormality in the outer packaging of a single bag, open the package, use the user to hold the end of the swab and cotton swab to collect the sample from the patient, gently rotate in the nasal cavity during the collection process, and do not use excessive force.

## Transportation and storage condition:

1. Stored in dry, clean and airy condition
2. Normal transportation


## Lot NO and expire date:

Lot NO and expire date, marked on the package

## Symbol description:



Manufacturer, Indicates the medical device manufacturer, as defined in Medical Device Regulation 2017/745 (EU)



Authorized representative in the European Community, Indicates the Authorized representative in the European Community



CE Mark: conforms to essential requirements of the Medical Device Regulation 2017/745 (EU)



Batch code, Indicates the manufacturer's batch code so that the batch or lot can be identified.



EO Sterilized



Keep in dry place , indicates a medical device need keep in dry place .



Do not use if package is damaged, indicates a medical device that should not be used if the package has been damaged or opened



Do not re-use, indicates a medical device that is intended for one use, or for use on a single patient during a single procedure



Caution, Indicates the need for the user to consult the instructions for use for important cautionary information, such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.




Use-by date, Indicates the date after which the medical device is not to be used



Name: Goodwood Medical Care Ltd

Address: 3-919 yongzheng street jinzhou district, dalian, China

Tel: 0411-82810067 Fax: 0411-82810468



Name: CMC Medical Devices & Drugs S.L.

Address :C/ Horacio Lengo No18, CP 29006, Málaga, Spain

Tel: +34951214054 Fax: +34952330100